



Mab to Hantavirus Seoul Serotype

Clone	R31
Category	Mouse monoclonal
Form/Purification	Hybridoma culture supernatant
Ig Class	IgG ₁
Immunogen	Nucleocapsid protein of Hantavirus R22 strain
Clinical Application/ Specificity	R31 reacts with an epitope present on the nucleocapsid protein of Hantavirus strain R22 (also called Seoul or SEO)
Application	<ul style="list-style-type: none">• ELISA• Immunofluorescence• Immunohistochemistry of frozen sections
Antigen Recognized in Species (tested so far)	Human, rat, mouse
Reconstitution Working Dilution	Reconstitute with 1 ml dist. water ready to use for immunohistochemistry
Storage	At 2-8°C
Quantity	1 ml (lyoph.); contains 0.09% NaN ₃ as preservative
Cat. No.:	R31-C